

OXYGEN

NIDEK NUVO NANO

PORTABLE OXYGEN CONCENTRATORS



Made In USA

Weighing at 2.1kg, the Nidek Nuvo Nano offers a peace of mind with up to 5 pulse oxygen flow setting and carry bag for portability. In-built smart safe audible and visual alerts prompt users and caregivers of any failure.



Quiet Operation



5 Setting Mode



Safe For Air Travels

Accessories



Nasal Cannula



Silicone Sleeve



Carry Bag



Power Adapter



2 Lithium-ion Battery



Replacement Intake Filter

Dimensions (L x W x H)	21.6 x 8.3 x 21.4 cm
Weight	2.2 kg
Electrical Requirements	100 - 240 VAC, 50-60 Hz
Oxygen Concentration	87% to 95.5%
Sound Level	< 40 dBA
Mode of Delivery	Pulse Mode Only
Oxygen Flow Volume	Setting 1: 210 ±15% ml/min Setting 2: 420 ±15% ml/min Setting 3: 630 ±15% ml/min Setting 4: 840 ±15% ml/min Setting 5: 1000 ±15% ml/min *±15% at Standard Temperature & Pressure
Battery Duration	Setting 1 : 3.5 hours Setting 2: 3.25 hours Setting 3: 2.5 hours Setting 4: 2.15 hours Setting 5: 1.5 hours
Operating Environment	Temperature Range: 5 - 40°C Altitude Range: Sea Level to 10,000ft (3,048m)
Storage Range	Temperature: 5°C to 35°C Humidity: Up to 90% (Non-Condensing)